

URETERAL ACCESS SHEATHS

NEGATIVE PRESSURE
& TIP FLEXIBLE



INNOVEX MEDICAL CO.,LTD.

2nd Floor, No.17 Building,
315 Qingda Rd, Pudong,
201201 Shanghai, P.R. China

Tel: +0086 021 50791830
Web: <http://www.innovexmed.com.cn> Zip: 201201
E-mail: global@innovexmed.com.cn




URETERAL ACCESS SHEATHS


NEGATIVE PRESSURE & TIP FLEXIBLE



- Control the negative pressure by adjusting the opening of vent



- Tip reinforcement to prevent deformation



Hydrophilic coating on sheath and dilator



- Sheath tip flexible

ORDER INFORMATION

| Product Code | Description | Specification |
|--------------|---|---------------|
| IVX-NS-1040 | tip flexible - negative pressure access sheath with scale, hydrophilic sheath and dilator | 10Fr x 40cm |
| IVX-NS-1050 | | 10Fr x 50cm |
| IVX-NS-1140 | | 11Fr x 40cm |
| IVX-NS-1150 | | 11Fr x 50cm |
| IVX-NS-1240 | | 12Fr x 40cm |
| IVX-NS-1250 | | 12Fr x 50cm |